



**Made In The USA!**

**Only until June 30th, 2017**

**Stock Up On Three Bio Plas Premium Lab Products At *20% Off!***

### G-Tube™ Premium Microcentrifuge Tubes:

Microcentrifuge Tubes without Compromise. We use only high quality medical grade ultra clear resins. G-Tubes are low-retention tubes allowing for maximum sample recovery. R-Nase & D-Nase Free, Non-pyrogenic endotoxin free. The G-Tube is built for high impact use and withstands over 30,000 x g when spun. The G-Tube will withstand temperature ranges of -80°C to 121 °C



### Hexa-Flex™ Safety Caps:

Fit a variety of tubes, including 10,12, 13, 16 and 18mm sizes of blood collecting and disposable culture tubes. The liquid-tight seal prevents sample evaporation. The Hexa-Flex™ Safety Cap is designed to hold under centrifugation, agitation, refrigeration and freezing. Lock securely onto, or into, the tube lip to prevent evaporation.



### Screw Cap Microcentrifuge Tubes:

Bio Plas Screw-Cap Micro Tubes are ideal for any sample processing or storage application and for working with valuable or hazardous cells and materials.

The extraordinary strength, safety, and reliability of the Bio Plas Tubes let you perform procedures that simply cannot be done with many other tubes. Spin it, store it, transport it (even by air). Boil it, irradiate it, immerse it in the vapor phase of liquid nitrogen. It will stand up to the test!

Bio Plas Screw-Cap Micro Tube's square threads won't strip and they are on the outside of the tube minimizing the risk of contaminating your sample or letting it contaminate you.



**BIO PLAS, INC**

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903  
USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com  www.bioplas.com