

New
Innovation

Drummond Portable Pipet-Aid® XL

**FREE T-Shirt
with Purchase!**
Choose from M, L or XL
Limited Supply

New Ergonomic Design Reduces Fatigue When Pipetting Under a Hood

- Longer lightweight handle lowers your arm position to reduce strain
- New quiet, more efficient pump extends operation period between charges
- Can be charged while in use to eliminate downtime
- Three-speed settings provide more precise control for both aspirate and dispense modes

Pipetting under a hood doesn't need to be a pain in the neck. The new Drummond Portable Pipet-Aid XL has a longer handle to reduce the amount of arm lift required to perform the same operation as with a conventional pipettor. The ergonomic design includes an adjustable sliding hand rest and a removable stand which enables the unit to be set down without contaminating the pipet. The Portable Pipet-Aid XL also features a new quieter, more powerful pump, and a new EZ Grip nosepiece. The unit is supplied complete with power supply charger, 4 extra filters, and a holster-type wall bracket.

DRUMMOND PORTABLE PIPET-AID® XL

Thomas No	Cat. No.	Description
1234X15	4-000-105	Portable Pipet-Aid XL 110V

Replacement Parts

Thomas No	Cat. No.	Description
7775L35	4-000-002	Rubber Insert (2/pk)
7775L75	4-000-051	TC Filters (5/pk)
0349P74	4-000-057	TC Gaskets (pair)
1233R82	4-000-036	Battery Pack (XL only)
1233R83	4-000-065	XL Power Supply/Charger, 110V
0349Q22	4-000-079	EZ Grip Nosepiece

Redemption Coupon

Offer Valid While Supplies Last

Complete the redemption coupon below and mail/fax/email it along with your Thomas Scientific proof of purchase to:

Thomas Scientific Marketing
1654 High Hill Road | Swedesboro, NJ 08085 USA
856.467.5167 (F)

Please send my FREE T-Shirt Size: M L XL (please circle one)

Name: _____

Company: _____

Street: _____

City: _____ State: _____

Zip/Postal: _____ Country: _____

Tel: _____

E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America, while supplies last. Limited to one redemptions per customer.



Compared to a conventional pipettor, the Portable Pipet-Aid XL lowers hand height required for the same task



The Developers of the Original Pipet-Aid



800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
www.thomassci.com • value@thomassci.com