

Safety SPOTLIGHT

Summer Edition

MICROFLEX
THE MOST TRUSTED NAME IN GLOVES

NeoPro® Green Chloroprene Gloves

- Unique Spring Green-colored gloves
- Textured fingertips for secure grip
- Latex- and powder-free to reduce skin irritation
- Improved durability and protection
- Soft feel and textured surface provides excellent comfort and tactile sensitivity
- Comfort and feel of natural latex
- Length is 9.5



Thomas No.	Description	Qty. Box	List Price	Sale Price	Qty. Case	List Price	Sale Price
5762E51	Neopro Green Gloves, X-Small	100	20.47	13.62	1000	184.20	136.20
5762E55	Neopro Green Gloves, Small	100	20.47	13.62	1000	184.20	136.20
5762E59	Neopro Green Gloves, Medium	100	20.47	13.62	1000	184.20	136.20
5762E63	Neopro Green Gloves, Large	100	20.47	13.62	1000	184.20	136.20
5762E67	Neopro Green Gloves, X-Large	100	20.47	13.62	1000	184.20	136.20

SafeGrip® Latex Gloves

- Polymer coated, powder-free, latex
- Blue, textured fingers
- Powder-free
- Ambidextrous
- Average thickness of no less than 12 mils
- Powder-free to help reduce the risk of skin irritation

Gloves have an extended, beaded cuff (290mm). Length: 11.4"; cuff thickness: 7.9 mils ±1.5; palm thickness: 11.0 mils ±0.8; finger thickness: 14.2 mils ±1.5.



Thomas No.	Description	Qty. Box	List Price	Sale Price	Qty. Case	List Price	Sale Price
5761W21	SafeGrip, Gloves, Small	50	23.16	15.08	500	208.40	150.80
5761W26	SafeGrip, Gloves, Medium	50	23.16	15.08	500	208.40	150.80
5761W31	SafeGrip, Gloves, Large	50	23.16	15.08	500	208.40	150.80
5761W36	SafeGrip, Gloves, X-Large	50	23.16	15.08	500	208.40	150.80

Thomas Scientific

We Believe You Are Important
— How Can We Help?



White Polypropylene Lab Coats with Elastic Wrists

- White
- Polypropylene
- No Pockets
- Elastic Wrists
- Sizes: SM - 5X



Thomas No.	Description	Qty. Case	List Price
1183A28	White Polypropylene Lab Coats with Elastic Wrists, S	30	33.27
1196L23	White Polypropylene Lab Coats with Elastic Wrists, M	30	35.12
1196L24	White Polypropylene Lab Coats with Elastic Wrists, L	30	35.12
1196L25	White Polypropylene Lab Coats with Elastic Wrists, XL	30	33.23
1196L26	White Polypropylene Lab Coats with Elastic Wrists, 2XL	30	33.79
1196L27	White Polypropylene Lab Coats with Elastic Wrists, 3XL	30	33.79
1196L28	White Polypropylene Lab Coats with Elastic Wrists, 4XL	30	33.79
1196L29	White Polypropylene Lab Coats with Elastic Wrists, 5XL	30	31.97

White Polypropylene Sleeves

- Polypropylene with elastic sleeves
- MATERIAL WEIGHT: 28 grams per square meter

ACCEPTANCES: All Safety Zone components comply with the provisions of the Federal Food, Drug, and Cosmetic Act which allow repeated use in direct contact with food.

Meets or Exceeds Class 100 Clean Room Standards



Thomas No.	Description	Qty. Case	List Price
1197U26	18" White Polypropylene Sleeve	200	32.01
1197U27	18" White Polyethylene Coated Sleeve	200	43.65

White Bouffant Caps

- Meets FFDC Requirements for Food Contact
- Light polypropylene
- Available in 19, 21 and 24 inch sizes
- Case of 1,000



Thomas No.	Description	Qty. Case	List Price
1197U17	19" White Bouffant Cap, Bags	1000	36.94
1197U18	21" White Bouffant Cap, Bags	1000	37.99
1197U19	24" White Bouffant Cap, Bags	1000	45.22



Think **Thomas** for All
of Your **Safety Needs!**

Order Now!

thomassci.com • global@thomassci.com
Inside US 800.524.0018 Outside US 856.467.2000
Fax 856.467.0512

Use Promo Code 0816 to receive the featured pricing. Valid until 12/31/16.



Think Thomas for All of Your Safety Needs!



New Pig

PIG® Absorbent Mat Pads

- Absorbs and retains oils, coolants, solvents and water
- Eight layers of 100% polypropylene are thermally bonded to make PIG Mat the strongest mat on the market; won't rip, tear or fray even when saturated
- Exclusive dimple pattern speeds wicking of liquid throughout mat for faster, easier cleanup
- Highly absorbent, fine-fiber construction won't leave liquids or fiber residue behind
- Absorbs and retains oils, coolants, solvents and water



Thomas No.	Description	Qty. Box	List Price	Qty. Case	List Price
1205U75	PIG® Absorbent Mat Pad, 20' x 15", Light-weight	—	—	200	103.38
1205U76	PIG® Absorbent Mat Pad, 20' x 15", Medium-weight	—	—	125	103.38
1205U77	PIG® Absorbent Mat Pad, 20' x 15", Heavy-weight	—	—	100	103.38
1186V32	PIG® Absorbent Mat Pad in Dispenser Box, 15" x 20", Heavyweight	100	106.31	—	—



Facemasks & Shields


With a variety of protective products and configurations, we have just the right one for your specific environment and/or application.



Thomas No.	Description	Qty. Case	List Price	Sale Price
1229E92	PFL N95 Particulate Respirator, Teal Stripe, Double Headband	210	206.16	173.33
1207K34	Facemask, Microbreathe™, earloops sonically welded on sides, white, cleanroom packaged (10 bags of 50)	500	177.29	160.46
1215J52	Facemask, AlphaGuard®, earloops sonically welded on sides, Blue, cleanroom packaged (10 bags of 50)	500	105.59	92.77
1215J54	Facemask, MVT, tie strings, sonically welded on sides, white, cleanroom packaged (10 bags of 50)	500	204.00	177.43
1215J53	Facemask, MVT, ear loops, sonically welded on sides, white, cleanroom packaged, (10 bags of 50)	500	137.80	127.08
1229E93	Full Face Shield w/ Comfort Band	100	263.22	217.33

Looking for signage?

Thomas Scientific offers a wide variety of signs from National Marker (NMC) to keep your workplace safe!



Custom signage also available, ask your Thomas Scientific Sales Representative.

Visit: www.thomassci.com, search: **National Marker**.



Astro Grip™ Powder-Free Nitrile Disposable Gloves




- 100 gloves per box
- Thickness: 6 mil
- Powder-free nitrile
- Dual sided scale grip
- Superior puncture & abrasion resistance

Thomas No.	Description	Qty. Box	List Price	Sale Price	Qty. Case	List Price	Sale Price
1188T56	Powder-Free Nitrile Gloves, M	100	15.49	12.50	1000	154.90	125.00
1188T57	Powder-Free Nitrile Gloves, L	100	15.49	12.50	1000	154.90	125.00
1188T58	Powder-Free Nitrile Gloves, XL	100	15.49	12.50	1000	154.90	125.00
1188Y66	Powder-Free Nitrile Gloves, 2XL	100	15.49	12.50	1000	154.90	125.00



Faucet Mounted Eyewash

Low-cost eyewash mounts easily to both standard and gooseneck faucets. For eyewash use, leave cold water supply on at all times. To operate unit pull actuator pin toward user. For faucet use, water flow is activated by the use of the diverter pin. Push the pin to the opposite side for water flow.



Thomas No.	Description	Qty. Each	List Price	Sale Price
1204M72	Faucet Mounted Eyewash	1	88.00	79.20

On-Site Portable Gravity-Fed Eyewash

- Holds 7 gallons of potable water
- Lower tank houses push-down eyewash tray
- SEI certified to ANSI Z358.1-2009



Thomas No.	Description	Qty. Each	List Price	Sale Price
1204M81	On-Site Portable Gravity-Fed Eyewash	1	351.00	315.90
1204M82	Safety Waste Water Cart for S19-921	1	415.00	373.50



Universal Spill Kit

Ampho-Mag™ is a rapid response powder that safely adsorbs and treats any spill and produces less heat and spattering than any other neutralizer available.



- Neutralizes acids and bases
- Color-coded
- Non-toxic
- Safe to handle
- Easy-to-store

Thomas No.	Description	Qty. Each	List Price	Sale Price
8238N05	Thomas® Universal Spill Kit	1	186.76	153.75



thomassci.com • global@thomassci.com

Inside US 800.524.0018
Outside US 856.467.2000
Fax 856.467.0512

Order Now!

Use Promo Code 0816
to receive featured pricing.
Valid through 12/31/16.

Think Thomas for All of Your Safety Needs!



First Aid Stations

First aid stations provide a large quantity and assortment of medical supplies for most minor injuries, aches and pains.

All stations are ANSI Z308, 1-2003 compliant



Thomas No.	Description	Qty. Each	List Price	Sale Price
1233T96	First Aid Station, 100 Person	1	184.55	153.78
1233T97	First Aid Station, 150 Person	1	265.51	221.25
1233T98	First Aid Station, 200 Person	1	431.71	381.68



Biodegradable Nitrile Disposable Gloves

- 100% nitrile with Eco Best Technology®
- 4 mil (0.10mm) thick
- 9.5 inches (240mm) long
- Zero natural rubber latex proteins
- Powder free
- Rolled cuff




Thomas No.	Description	Qty. Box	List Price	Sale Price	Qty. Case	List Price	Sale Price
1182Y72	Biodegradable Nitrile Glove, Green, Powder-Free, XS	100	13.20	12.00	1000	132.00	120.00
1182Y73	Biodegradable Nitrile Glove, Green, Powder-Free, S	100	13.20	12.00	1000	132.00	120.00
1182Y74	Biodegradable Nitrile Glove, Green, Powder-Free, M	100	13.20	12.00	1000	132.00	120.00
1182Y75	Biodegradable Nitrile Glove, Green, Powder-Free, L	100	13.20	12.00	1000	132.00	120.00
1182Y76	Biodegradable Nitrile Glove, Green, Powder-Free, XL	100	13.20	12.00	1000	132.00	120.00
1182Y77	Biodegradable Nitrile Glove, Green, Powder-Free, XXL	100	13.20	12.00	1000	132.00	120.00



Multi-Purpose Lab Soakers

Spilfyter® Soakers provide absorption protection against spills and cushioning to minimize breakage.

Soakers are made with multiple layers of 100% recycled tissue and a plastic poly backing for diaper-like characteristics. Heavy duty models feature durable fabric top.



Thomas No.	Description	Qty. Case	List Price	Sale Price
1233K17	Lab Soaker, Light Duty Roll, 400' X 20"	2	180.17	165.22
1233K18	Lab Soaker, Light Duty Pad, 20" X 16'	250	82.70	75.88
1233K19	Lab Soaker, Heavy Duty Roll w / Fabric Top, 200' X 20"	2	145.96	133.82
1233K20	Lab Soaker, Heavy Duty Pad w / Fabric Top, 20" X 16'	200	108.01	99.12



Auto-Retractable Box Cutter

- Auto-retractable for extra safety
- High visibility green rubberized interior grip
- Includes one double-sided ceramic blade
- Stays sharp 10X longer than steel



Thomas No.	Description	Qty. Case	List Price	Sale Price
1191X86	Box Cutter, Auto-Retractable, Ceramic	6	195.95	152.62

Safety Cutter

- Perfect for "knife free" facilities
- Clean-room friendly no oily blades!
- Ceramic blade lasts longer than steel
- Stays sharp and never rusts



Thomas No.	Description	Qty. Case	List Price	Sale Price
1191X83	Safety Cutter, Ceramic, Carded	12	117.51	91.42
1191X84	Safety Cutter, Ceramic, Pillow Pack	48	470.03	365.67



Simple Green Safety Towels®

- Specially designed with the industrial cleaning strength of Crystal Simple Green infused into thermally bonded, dual textured, melt-blown, lint-free polypropylene towels
- Cleaning towels are Pre-moistened with non-toxic and biodegradable Crystal Simple Green®
- Excellent wicking/fluid-holding capacity, are highly resistant to rot and mildew and have an abrasive side for tough scrubbing.
- Multi-purpose cleaning towels can handle just about any clean-up challenge you may face — and do so with the highest regard for human safety and environmental protection



Quick, effective and versatile Simple Green Safety Towels are the solution for your toughest, multi-purpose cleaning tasks.

Thomas No.	Description	Qty. Case	List Price	Sale Price
1200X47	Simple Green Safety Towels, 75 Count	6	154.85	120.51

Lemon Scent Simple Green® All-Purpose Cleaner

- Original non-toxic, biodegradable, concentrated formula
 - Available in lemon scent
- The concentrated formula can be diluted with water to remove grease, soap scum, coffee and juice stains, ink, lipstick and more from any washable, indoor or outdoor surface. Lemon Scent is ideal for kitchens and baths!



Thomas No.	Description	Qty. Case	List Price	Sale Price
1200X70	Lemon Scent Simple Green All-Purpose Cleaner, 24 oz	12	132.68	111.38
1200X71	Lemon Scent Simple Green All-Purpose Cleaner, 1 gallon	6	136.75	114.80



thomassci.com • global@thomassci.com

Inside US 800.524.0018
Outside US 856.467.2000
Fax 856.467.0512

Order Now!

Use Promo Code 0816
to receive featured pricing.
Valid through 12/31/16.

Order Now!

thomassci.com • global@thomassci.com

Inside US 800.524.0018
Outside US 856.467.2000
Fax 856.467.0512

Use Promo Code 0816
to receive featured pricing.

Valid through 12/31/16.



Bouton® Optical The Scout Visitor Specs

- Polycarbonate unilens
- Over-the-glass design fits over all prescription eyewear
- Can also be worn by itself as a regular safety glass
- Molded nose bridge design to fit the vast majority of users
- Indirect vent temples to reduce fogging
- 4 Base curve lens



Thomas No.	Description	Qty. Pack	List Price	Sale Price	Qty. Case	List Price	Sale Price
1200W90	The Scout Eyewear, Visitor Specs, Clear Temple, Clear Lens, Anti-Scratch	36	70.71	47.78	144	226.28	191.12
1200W91	The Scout Eyewear, Visitor Specs, Gray Temple, Gray Lens, Anti-Scratch	36	70.71	50.38	144	226.28	201.52
1200W92	The Scout Eyewear, Visitor Specs, Clear Temple, Clear Lens, Uncoated	36	55.91	42.55	144	178.92	170.20



Sure-Grip® EX Flammable Safety Cabinets

- Protect workers, reduce fire risks, and improve productivity by storing flammable fuels and chemicals in code-compliant safety cabinets
- High-performance, self-latching doors close easily and securely for maximum protection under fire conditions
- Sure-Grip® EX cabinets are constructed of sturdy 18-gauge (1-mm) thick double-wall, welded steel with 1-1/2" (38-mm) of insulating air space for fire resistance
- For added security, U-Loc™ handle comes with a cylinder lock and key set or accepts an optional padlock



Thomas No.	Description	Qty. Each	List Price	Sale Price
0501A44	Sure-Grip EX Flammable Safety Cabinet, Yellow, 30 gal., 1 shelf, 1 bi-fold self-close door	1	1482.97	1235.82

Sure-Grip® EX Slimline Flammable Safety Cabinets

- Use alone or alongside larger cabinets for expanded storage when floor space is limited
- For added security, U-Loc™ handle comes with a cylinder lock and key set or accepts an optional padlock
- Patented SpillSlope® galvanized steel shelves direct spills to back and bottom of a leakproof 2" (51-mm) bottom sump



Thomas No.	Description	Qty. Each	List Price	Sale Price
1185Q96	Sure-Grip EX Slimline Flammable Safety Cabinet, Yellow, Cap. 22 gallons, 3 shelves, 1 manual-close door	1	927.32	772.77
1184F86	Sure-Grip EX Slimline Flammable Safety Cabinet, Yellow, Cap. 22 gallons, 3 shelves, 1 self-closing door	1	978.03	815.03



GHS Compliant, Multi-Lingual Wash Bottles

- GHS Compliant, Pre-Printed Wash Bottles
- Safety-Vented (1190V90 - V97)
- Non-Vented (1190V98 - V99)



Thomas No.	Description	Qty. Box	List Price	Sale Price
1190V90	Wash Bottle, Acetone, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Red (Pack/5)	5	39.40	36.37
1190V91	Wash Bottle, Acetone, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Red (Each)	1	10.20	9.00
1190V92	Wash Bottle, Methanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Green (Pack/5)	5	39.40	36.37
1190V93	Wash Bottle, Methanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Green (Each)	1	10.20	9.00
1190V94	Wash Bottle, Isopropanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Yellow (Pack/5)	5	39.40	36.37
1190V95	Wash Bottle, Isopropanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Yellow (Each)	1	10.20	9.00
1190V96	Wash Bottle, Ethanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Orange (Pack/5)	5	39.40	36.37
1190V97	Wash Bottle, Ethanol, 500mL, LDPE, Multi-Lingual, Safety Vented, GHS Compliant, Orange (Each)	1	10.20	9.00
1190V98	Wash Bottle, Distilled Water, 500mL, LDPE, Multi-Lingual, GHS Compliant, White (Pack/5)	5	30.80	28.43
1190V99	Wash Bottle, Distilled Water, 500mL, LDPE, Multi-Lingual, GHS Compliant, White (Each)	1	8.50	7.50



SemperShield® Nitrile Gloves

- Extra thickness and an extended cuff
- Excellent chemical resistance
- Tested for use with chemotherapy drugs
- Manufactured without natural rubber
- Powder-Free



Thomas No.	Description	Qty. Box	List Price	Qty. Case	List Price
1186P92	SemperShield Nitrile Glove, 5 mil, Powder-Free, 12" Extended Cuff, S	50	6.50	500	65.00
1186P93	SemperShield Nitrile Glove, 5 mil, Powder-Free, 12" Extended Cuff, M	50	6.50	500	65.00
1186P94	SemperShield Nitrile Glove, 5 mil, Powder-Free, 12" Extended Cuff, L	50	6.50	500	65.00
1186P95	SemperShield Nitrile Glove, 5 mil, Powder-Free, 12" Extended Cuff, XL	50	6.50	500	65.00
1186P96	SemperShield Nitrile Glove, 5 mil, Powder-Free, 12" Extended Cuff, XXL	50	6.50	500	65.00



We Believe You Are Important
— How Can We Help?

Inside US 800.524.0018
Outside US 856.467.2000 • Fax 856.467.0512
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • global@thomassci.com

