

Receive FREE
Boekel Products!



Purchase a Shaking Water Bath

Thomas No. 1229U55

Get a universal flask tray
Free!

- Switch from orbital to reciprocating mode by simply reversing the tray
- 3 Point factory programmed temperature calibration throughout temperature range
- PID Controller for easy temperature selection, rapid heat up and excellent stability + 0.1°C
- Programmable timer up to 99 hrs, 59 min – Bath stops heating when timed operation is completed
- Adjustable speed control from 20-200 RPM (in 1 RPM increments)
- Convenient drain port with quick disconnect allows for easy cleaning of the bath

Purchase a Mini Tube Rotator

Comes with one carousel Thomas No. 1217H25

Get remaining 2 carousels *Free!*

- Variable Speed: 4-18 rpm
- LCD display
- Includes Thomas No. 1217H26 Carousel for thirty-six 1.5/2.0ml tubes
- Provides complete 360° Rotation
- An adjustable 0-90° Mixing Angle



Please see page 2 for redemption form



We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • value@thomassci.com





Purchase a **Shaking Water Bath**

Get a universal flask tray *Free!*

Thomas No. 1229U57

Purchase a **Mini Tube Rotator**

Get remaining 2 carousels *Free!*

Please select your FREE carousels from list below.

- 1217H26
- 1217H27
- 1217H28

Offers subject to change without notice. Please allow 4-6 weeks for delivery of free products when applicable. Promotion expires December 31, 2017. All redemption forms must be received no later than January 15, 2018 to be valid.

Please complete the account information below, and fax or email, along with Thomas Scientific Proof-of-Purchase, to:
Thomas Marketing: Fax 856-467-5167 Email Marketing@thomassci.com

Account Information

Name _____ Date _____

Organization _____

Account No. _____

Ship To Address _____

Room/Bldg./Dept. _____

City _____ State _____ Zip _____

Phone _____

Email _____



We Believe You Are Important
— How Can We Help?

800.345.2100 • 800.345.5232 fax
PO Box 99 • Swedesboro, NJ 08085
thomassci.com • value@thomassci.com

