

Multi-Tube Vortexer

Hands-free vortexing of up to 50 tubes

- Automatic timer
- Electronic speed control

Thomas multi-tube vortexer can handle up to 50 tubes. Speed is adjustable from 1200 to 2400 rpm with electronic speed control. Built-in timer is adjustable from one second to one minute or continuous operation. Various tube sizes can be mixed by simply changing the foam tube racks. Unit easily adjusts to different tube heights. Solid steel posts secure tube support for safe operation. Heavy-duty base features a scratch-resistant enamel finish.



Dimensions: 17.25" W x 9" D x 16" H
Shipping weight: 45 lbs.

Thomas No.	Description	Each
8294C01	Multi-Tube Vortexer, 120V, 50/60 Hz	1
8294C03	Multi-Tube Vortexer, 230V, 50/60 Hz	1
8294C11	10 mm Tube Rack, Foam, Gray	1
8294C13	12 mm Tube Rack, Foam, Blue	1
8294C15	13 mm Tube Rack, Foam, Yellow	1
8294C17	16 mm Tube Rack, Foam, Green	1
8294C19	Post Extension Kit	1