

# NEST Biotechnology Co.,Ltd

## Certificate of Analysis of Disposable Sampler

Product Name	Disposable sampler	Product No.	202001	Lot No.	030120NSM02
No.	Test items	Standards			Results
1	Appearance	No maculas, bubble and scratch, etc			Pass
2	Net Content	$\geq 2\text{ml}$			Pass
3	PH	7.2-7.8			Pass
4	The Effect of Preservation	The variation coefficients of CT values were all less than 5%			Pass
5	Dimensions	Meet the requirement of the drawing			Pass
6	Seal Test	Without disclosure of G12 in 3 minutes under -30kpa~-29Kpa vacuum			Pass
7	Vibration Test	No disclosure after circular shaker 210rpm, 2hours			Pass
8	Freezing Test	No cracking/breakage/leakage/deformation at -80℃			Pass
9	Effect of Screw Cap	The tube body and cap should not be damaged, and the cap should be properly tightened			Pass
10	Packaging	Packaging materials and quantity is correct, packaging is intact			Pass
Conclusion		OK			
Note:"Results":OK or NG					
Address: NO.530,XIDA ROAD,NEW DISTRICT,WUXI		Drawn:	赵静	Approved:	杨洋
TEL: (+86) 0510-88550100					
FAX: (+86) 0510-88550105		Date:	2020.3.2	Date:	2020.3.2
http://www.cell-nest.com					

**NEST**

F-BD-032

Edition:02

