

# NON-CONTACT

## INFRARED FOREHEAD THERMOMETERS

*FDA CLEARED*



Non-Contact Infrared Forehead Thermometer

Power Supply: 2 AAA Batteries Included

Celsius/Fahrenheit Switchable

Range: 89.6/109.2F & 35.0/42.9C

Accuracy: +/- 0.4F & 0.2C

Resolution: 0.1F & 0.1C

Approximate Size of Product: 95mm x 60mm x 180mm

Approximate Weight of Product: 146g

Memory: 32 Measured values

Automatic Shutdown when not in use

*PART NUMBER: GP-300*

*Take Adult & Child Forehead  
Temperatures Safely & Quickly with  
No-Touch Infrared Technology.*

*We supply Forehead Thermometers  
with Proven Accuracy.*

**KESSLER**  
Health Products

(631) 841-5500  
info@kesslerusa.com  
www.kesslerusa.com