

10" x 12" Clear Face Shields

Affordable protection for healthcare, sanitation, industrial, and food service professionals.



Face shields provide frontline personnel basic protection from sprayed liquid particles and can be wiped clean for safe reuse. Our face shields are constructed with a latex-free Velcro strap for the most comfortable fit anywhere.

Clear laminate made of lightweight .0052" thickness OPE optical polyester

- Polyester laminate dimensions: 10" tall x 12" wide, with comfort cut corners
- Laminate extends 0.5" above foam head strip and 8.5" below headstrip
- Fog-resistant
- Crease-resistant

1" x 1.5" thick foam head strip (strip is 9" long)

- Easily accommodates eyeglasses and protective goggles
- Compatible with stethoscope use

Adjustable Velcro headstrap

- Velcro offers the most comfortable fit of any strap. Tighten in windy outdoor conditions, or loosen indoors
- Completely adjustable, unlike elastic bands
- Does not pinch like "eyeglass" type plastic frames
- One size fits all

Additional features:

- Latex-free
- Can be cleaned with bleach, ammonia, or peroxide wipes without streaking or discoloring
- Can be cut to any custom size with common household scissors
- Non-sterile. Recommended to be used in combination with additional protective gear according to health guidelines.
- 144-count carton. Minimum order one carton. Bulk packed in one plastic bag.
- Carton weight 12 lbs.
- Pallet size: 25 cartons

Thomas No.	Description	Quantity	Price
20A00G602	Face Shield, Clear 10" x 12" Single-Use	144/Case	\$1,058.50

Disclaimer: Pricing and offers subject to change without notice.